



For braces placed on the inside of the teeth.

TePe Universal Care™

This brush has an angled neck which makes it easy to reach braces on the inside.

TePe PlaqSearch™

Use a disclosing tablet which colours plaque to check if you have brushed properly. PlaqSearch shows both new and old plaque. Ask your dentist or dental hygienist before use.



We care for healthy smiles

Interdental brushes

Original	Extra soft	Angle
0.4 mm <input type="checkbox"/>	0.45 mm <input type="checkbox"/>	0.4 mm <input type="checkbox"/>
0.45 mm <input type="checkbox"/>	0.5 mm <input type="checkbox"/>	0.45 mm <input type="checkbox"/>
0.5 mm <input type="checkbox"/>	0.6 mm <input type="checkbox"/>	0.5 mm <input type="checkbox"/>
0.6 mm <input type="checkbox"/>	0.7 mm <input type="checkbox"/>	0.6 mm <input type="checkbox"/>
0.7 mm <input type="checkbox"/>	0.8 mm <input type="checkbox"/>	0.7 mm <input type="checkbox"/>
0.8 mm <input type="checkbox"/>	1.1 mm <input type="checkbox"/>	0.8 mm <input type="checkbox"/>
1.1 mm <input type="checkbox"/>		
1.3 mm <input type="checkbox"/>		
1.5 mm <input type="checkbox"/>		

© TePe Mundhygienprodukt AB, 2019 / 173INT.

TePe EasyPick™

- Orange XS/S
- Turkos M/L

Toothbrushes

- Select Compact, x-soft
- Select Compact, soft
- Select, x-soft
- Select, soft
- Supreme

Special brushes

- Orthodontic Brush
- Interspace, soft
- Interspace, medium
- Compact Tuft
- Universal Care

Additional products

- PlaqSearch

BR171069INT/vers4

Made in Sweden

Tand Co., Ltd.

7 Soi On Nut 1/1, Sukhumvit 77
Sukhumvit Road, Watthana
Bangkok 10110
Tel +66 2 092 2031-2
info@tandthailand.com
www.tandthailand.com



Caring for your braces



Tips for a healthy smile

- Brush with a soft toothbrush and fluoride toothpaste twice a day.
- Use a special brush for additional cleaning in difficult to reach areas.
- Avoid sweet drinks and snacks between meals.
- For extra caries protection rinse with fluoride after toothbrushing.

Take care of your smile

Good oral hygiene is the key to a healthy mouth throughout life. With braces tooth cleaning is more complicated as plaque quickly builds up around brackets, bands and wires. Daily cleaning with a regular toothbrush is important to prevent cavities, gum inflammation and bad breath.

TePe Select™ Compact

Select Compact has a small brush head for better access in difficult to reach areas.

TePe Select™

Select has a slightly bigger brush head than Select Compact.

TePe Supreme™



This toothbrush has filaments in two levels. The longer filaments easily clean under the wire.

TePe Orthodontic Brush

The extra narrow brush head enables cleaning close to brackets and under the wire.



Brushing technique




You may use an additional product for areas that are difficult to reach:




TePe Interspace™

This angled brush with a slim and pointed tuft is ideal for cleaning critical areas close to the brackets and under the wire. Interchangeable brush heads.



TePe Compact Tuft™

The angled brush with a small, rounded tuft is suitable for cleaning around the brackets and along the gum line.



Under the wire.

Between the teeth.

TePe Interdental Brushes

Interdental brushes are useful for cleaning under the wire and between the teeth, if there is enough space.

TePe EasyPick™

Easy cleaning between the teeth under the wire or retainer.

